



3 professional tests for your practice




3 rapid tests based on finger tip blood

Biocard™ Helicobakter Pylori Test, CeliacTest and CanDia5® are scientifically proven and reliable.


All tests are compared to accepted reference methodologies and perform at the level of other pathology tests or better.

**CanDia5®
Candida Test**




| | |
|-------------|-----|
| Specificity | 90% |
| Sensitivity | 78% |

**Biocard™
Celiac Test**



| | |
|-------------|-------|
| Specificity | 93.5% |
| Sensitivity | 96.7% |

**Biocard™
Helicobacter
Pylori Test**



| | |
|-------------|-------|
| Specificity | 97.8% |
| Sensitivity | 93.2% |

For more information on these test please contact

CanDia5[®]

Candida Rapid Test




CanDia5[®] has been designed specifically for practitioners enabling them to offer laboratory accuracy for the detection of Candida at the time of patient consultation.

CanDia5[®] offers:

- **Rapidity** – the result is available in 5 minutes
- **Reliability** – CanDia5[®] has been rigorously tested for accuracy and reproducibility
- **Convenience** – the test procedure requires only 5 easy steps

The **CanDia5[®] Candida Rapid Test** is an immunochromatographic test kit that accurately detects circulating antibodies specific to a current or recent Candida infection. With only one drop of fingertip blood, this self-contained test kit detects the four most common species of Candida antibodies.



Biocard[™] Celiac Test

For the rapid detection of Coeliac Disease

1 to 2 in 100 individuals suffer from Coeliac Disease (1–2 %).
Only 10 to 20% of CD sufferers are diagnosed.

Coeliac Disease can have serious long term medical implications such as:

- | | | |
|------------------|-------------------|-------------------|
| ■ Stunted growth | ■ Chronic fatigue | ■ Anaemia |
| ■ Osteoporosis | ■ Infertility | ■ Lymphoma |
| ■ Liver disease | ■ Depression | ■ Immune weakness |

Approximately 250,000 Australians suffer from Coeliac Disease, 200,000 to 225,000 remain undiagnosed!
 CD is the most under-diagnosed disease in Australia.

The **Biocard[™] Celiac Test** is fast (5mins), easy to use and as accurate as the best pathology tests. The BCT is perfect for screening for CD as well as for monitoring diet compliance of CD sufferers.



Biocard[™] Helicobacter Pylori Test

For the rapid detection of H Pylori

Approximately 50% of the population is infected with Helicobacter Pylori.

H Pylori can cause gastritis in adults and children plus it causes ulcers. Persons infected with H Pylori have a 2 to 6 fold increased risk of developing gastric cancer and lymphoma as compared to uninfected persons.

NEW!

The **Biocard[™] Helicobacter Pylori Test** is a fingertip blood test, delivers a result in 5 minutes and offers an excellent accuracy of 96.2%

